

Product Datasheet

Cd2 Antibody (orb1272061)

Description	Cd2 Antibody
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC
Target	Cd2
Preservatives	Phosphate-buffered aqueous solution, pH 7.2.
Form/Appearance	liquid
Concentration	batch dependent
Storage	The product should be stored undiluted at 4°C . Do not freeze.
Note	For research use only
Isotype	Rat IgG2b, lambda
Clonality	Monoclonal
Uniprot ID	P08920
Expiration Date	12 months from date of receipt.