


Product Datasheet

CD2 Antibody (orb1272062)

Description	CD2 Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FA, FC
Target	CD2
Preservatives	Phosphate-buffered aqueous solution, pH 7.2.
Form/Appearance	liquid
Concentration	batch dependent
Storage	The product should be stored undiluted at 4°C . Do not freeze.
Note	For research use only
Isotype	Mouse IgG1, kappa
Clonality	Monoclonal
Uniprot ID	P06729
Expiration Date	12 months from date of receipt.



Human
peripheral
blood
lymphocytes
were ...