

Product Datasheet

Cd81 Antibody (orb1272064)

Description	Cd81 Antibody
Species/Host	Hamster
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC, IP, WB
Target	Cd81
Preservatives	Phosphate-buffered aqueous solution, pH 7.2.
Form/Appearance	liquid
Concentration	batch dependent
Storage	The product should be stored undiluted at 4°C . Do not freeze.
Note	For research use only
Isotype	Armenian Hamster IgG
Clonality	Monoclonal
Uniprot ID	P35762
Expiration Date	12 months from date of receipt.