

Product Datasheet

Itgb3 Antibody (orb1272065)

Description	Itgb3 Antibody
Species/Host	Hamster
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FA, FC
Target	Itgb3
Preservatives	Phosphate-buffered aqueous solution, pH 7.2.
Form/Appearance	liquid
Concentration	batch dependent
Storage	The product should be stored undiluted at 4°C . Do not freeze.
Note	For research use only
Isotype	Armenian Hamster IgG
Clonality	Monoclonal
Uniprot ID	O54890
Expiration Date	12 months from date of receipt.