

Product Datasheet

ITGB3 Antibody (FITC) (orb1272069)

Description	ITGB3 Antibody (FITC)
Species/Host	Mouse
Reactivity	Human
Conjugation	FITC
Tested Applications	FC
Target	ITGB3
Preservatives	Phosphate-buffered aqueous solution, ≤0.09% Sodium azide, may contain carrier protein/stabilizer, pH 7.2.
Form/Appearance	liquid
Concentration	batch dependent
Storage	The product should be stored undiluted at 4°C and should be protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	mouse IgG1, kappa
Clonality	Monoclonal
Uniprot ID	P05106
Expiration Date	12 months from date of receipt.