

Product Datasheet

MSK2 (phospho-Thr568) antibody (orb127207)

Description	Rabbit Polyclonal to MSK2 (phospho-Thr568).
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	O75676
Dilution Range	WB: 1:500-1:3000, ELISA: 1:10000
Expiration Date	12 months from date of receipt.