

Product Datasheet

ITGB3 Antibody (orb1272071)

Description	ITGB3 Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FA, FC, IHC, WB
Target	ITGB3
Preservatives	Phosphate-buffered aqueous solution, ≤0.09% Sodium azide, may contain carrier protein/stabilizer, pH 7.2.
Form/Appearance	liquid
Concentration	batch dependent
Storage	The product should be stored undiluted at 4°C . Do not freeze.
Note	For research use only
Isotype	mouse IgG1, kappa
Clonality	Monoclonal
Uniprot ID	P05106
Expiration Date	12 months from date of receipt.