

Product Datasheet

p70 S6 Kinase beta (phospho-Ser423) antibody (orb127209)

Description Rabbit Polyclonal to p70 S6 Kinase beta (phospho-Ser423).

Species/Host Rabbit

Reactivity Human, Mouse

Conjugation Unconjugated

Tested Applications WB

Immunogen C-terminal

Concentration 1mg/ml


Note For research use only

Clonality Polyclonal

Uniprot ID **Q9UBS0**

Dilution Range WB: 1:500-1:3000, ELISA: 1:1000

Expiration Date 12 months from date of receipt.



Western
blot
analysis of
extracts
from K...