


Product Datasheet

ITGAM Antibody (PE/Cy5) (orb1272118)

Description	ITGAM Antibody (PE/Cy5)
Species/Host	Rat
Reactivity	Human, Mouse
Conjugation	PE/Cy5
Tested Applications	FC
Target	ITGAM
Preservatives	Phosphate-buffered aqueous solution, ≤0.09% Sodium azide, may contain carrier protein/stabilizer, pH 7.2.
Form/Appearance	liquid
Concentration	batch dependent
Storage	The product should be stored undiluted at 4°C and should be protected from prolonged exposure to light. Do not freeze.
Note	For research use only
Isotype	Rat IgG2b, kappa
Clonality	Monoclonal
Uniprot ID	P11215
Expiration Date	12 months from date of receipt.



C57Bl/6
bone
marrow cells
were stained
w...