


## Product Datasheet

### SHB antibody (orb127226)

|                            |                                 |
|----------------------------|---------------------------------|
| <b>Description</b>         | Rabbit Polyclonal to SHB.       |
| <b>Species/Host</b>        | Rabbit                          |
| <b>Reactivity</b>          | Human, Mouse                    |
| <b>Conjugation</b>         | Unconjugated                    |
| <b>Tested Applications</b> | WB                              |
| <b>Immunogen</b>           | C-terminal                      |
| <b>Concentration</b>       | 1mg/ml                          |
| <b>Note</b>                | For research use only           |
| <b>Clonality</b>           | Polyclonal                      |
| <b>Uniprot ID</b>          | <a href="#">Q15464</a>          |
| <b>Dilution Range</b>      | WB: 1:500-1:3000, ELISA: 1:5000 |
| <b>Expiration Date</b>     | 12 months from date of receipt. |



Western blot analysis of HuvEc cells using SHB antibody.