


Product Datasheet

CCL20 Antibody (orb1272363)

Description

CCL20 Antibody

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Produced from sera of rabbits pre-immunized with highly pure (>98%) recombinant hMIP-3-alpha (human Macrophage Inflammatory Protein-3-alpha).
Target	CCL20
Form/Appearance	Lyophilized
Concentration	batch dependent
Storage	MIP-3a antibody is stable for at least 2 years from date of receipt at -20°C. The reconstituted antibody is stable for at least two weeks at 2-8°C. Frozen aliquots are stable for at least 6 months when stored at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P78556
NCBI	P78556
Dilution Range	ELISA:To detect hMIP-3-alpha by direct ELISA (using 100 µL/well antibody solution) a concentration of at least 0.5 µg/mL of this antibody is required. This antigen affinity purified antibody, in conjunction with compatible secondary reagents, allows the detection of 0.2 - 0.4 ng/well of recombinant hMIP-3-alpha.Sandwich:To detect hMIP-3α by sandwich ELISA (using 100 µL/well antibody solution) a concentration of 0.5 - 2.0 µg/mL of this antibody is required. This antigen affinity purified antibody, in conjunction with our Biotinylated Anti-Human MIP-3αas a detection antibody, allows the detection of at least 0.2 - 0.4 ng/well of recombinant hMIP-3α. Western Blot:To detect hMIP-3α by Western Blot analysis this antibody can be used at a concentration of 0.1-0.2 µg/mL. Used in conjunction with compatible secondary reagents the detection limit for recombinant hMIP-3α is 1.5-3.0 ng/lane, under either reducing or non-reducing conditions.
Expiration Date	12 months from date of receipt.



To detect Human MIP-3-alpha by sandwich ...



To detect HumanMIP-3-alpha by Western Bl...



To detect HumanMIP-3-alpha by Western Bl...