



Product Datasheet

PTPRA antibody (orb12726)


Description	Rabbit polyclonal antibody to PTPRA
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human PTPRA(NP_543031.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	120kDa
Uniprot ID	P18433
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human stomach tissue.



Immunohistochemical staining of human kidney tissue.



Immunohistochemical staining of human colon tissue.