


Product Datasheet

IL17B Antibody (orb1272618)


Description

IL17B Antibody


Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	E.coli derived Recombinant Human IL-17B
Target	IL17B
Form/Appearance	Lyophilized
Concentration	batch dependent
Storage	The lyophilized antibody is stable for at least 2 years from date of receipt at -20°C. The reconstituted antibody is stable for at least two weeks at 2-80°C. Frozen aliquots are stable for at least 6 months when stored at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	ELISA: Sandwich: To detect hIL-17B by sandwich ELISA (using 100 µL/well antibody solution) a concentration of 0.5 - 2.0 µg/mL of this antibody is required. This antigen affinity purified antibody, in conjunction with Biorbyt's Polyclonal Biotinylated Anti-Human IL-17B (38-173) as a detection antibody, allows the detection of at least 0.2 - 0.4 ng/well of recombinant hIL-17B. Western Blot: To detect hIL-17B by Western Blot analysis this antibody can be used at a concentration of 0.1 - 0.2 µg/mL. Used in conjunction with compatible secondary reagents the detection limit for recombinant hIL-17B is 1.5 - 3.0 ng/lane, under either reducing or non-reducing conditions.
Clonality	Polyclonal
Uniprot ID	Q9UHF5
NCBI	Q9UHF5
Expiration Date	12 months from date of receipt.



To detect hIL-17B by sandwich ELISA (usi...



To detect hIL-17B by Western Blot analys...



To detect hIL-17B by Western Blot analys...