

Product Datasheet

SARS-CoV-2 (COVID-19) NSP1 polyclonal antibody (orb1272805)

Description SARS-CoV-2 (COVID-19) NSP1 polyclonal antibody

Species/Host Rabbit

Reactivity Virus

Conjugation Unconjugated

Tested Applications ELISA

Immunogen Recombinant SARS-CoV-2 NSP1(Met1-Gly180)

Preservatives PBS, pH 7.4, containing 0.05% proclin300, 50% glycerol.

Form/Appearance Liquid

Concentration batch dependent

Storage Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8 °C for one week. Store at -20 to -80 °C for twelve months from the date of receipt.

Note For research use only


Application notes Elisa:1:4000~1:8000

Isotype IgG

Clonality Polyclonal

Dilution Range Elisa:1:4000~1:8000

Expiration Date 12 months from date of receipt.



Various lysates were subjected to SDS-PAGE