

Product Datasheet

IL-33 Antibody [IL33305B] (orb1272817)

Description

IL-33 Antibody [IL33305B]

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, FA, IHC, IP, WB

Applications

Immunogen

Recombinant human IL-33.

Preservatives

Liquid. 0.2µm-filtered solution in PBS, pH 7.4.
Contains no preservatives.

Form/Appearance

Liquid

Concentration

1 mg/ml

Storage

Stable for at least 1 year after receipt when stored
at -20°C.

Note

For research use only

Application notes

ELISA: (direct or indirect: 1:2'000-1:5'000),
Functional Application (Blocking),
Immunohistochemistry (paraffin sections (1:50-
1:200)), Immunoprecipitation Western blot:
(1:2'000-1:5'000 using ECL. Suggested blocking
and dilution buffer is PBST with 0.05% Tween 20
and 5% skim milk. Suggested incubation time is 1
hour at room temperature). Optimal conditions
should be determined individually for each
application.

Isotype

IgG2ak

Clonality


Monoclonal

Dilution Range


ELISA: (direct or indirect: 1:2'000-1:5'000),
Functional Application (Blocking),
Immunohistochemistry (paraffin sections (1:50-
1:200)), Immunoprecipitation Western blot:
(1:2'000-1:5'000 using ECL. Suggested blocking
and dilution buffer is PBST with 0.05% Tween 20
and 5% skim milk. Suggested incubation time is 1
hour at room temperature). Optimal conditions
should be determined individually for each
application.

Expiration Date


12 months from date of receipt.



Immunohistochemical
staining of human ti...



Inhibition of binding
of recombinant
hum...



Immunoprecipitation
of recombinant hIL-
3...