

Product Datasheet

ANGPTL3 Antibody (orb1272858)

Description

ANGPTL3 Antibody

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant human ANGPTL3.
Target	ANGPTL3
Preservatives	Liquid. 0.2um-filtered solution in PBS, pH 7.4. Contains no preservatives.
Form/Appearance	Liquid
Concentration	1 mg/ml
Storage	Stable for at least 6 months after receipt when stored at -20°C.
Note	For research use only
Application notes	ELISA: (direct and indirect: 1:2,000-1:5,000). Western Blot: (1:2,000-1:5,000 using ECL; suggested blocking and dilution buffer is PBST containing 0.05% Tween 20 and 5% skim milk; suggested incubation time is 1 hour at room temperature). Optimal conditions must be determined individually for each application.
Clonality	Polyclonal
Uniprot ID	A0JLS0
NCBI	A0JLS0
Dilution Range	ELISA: (direct and indirect: 1:2,000-1:5,000). Western Blot: (1:2,000-1:5,000 using ECL; suggested blocking and dilution buffer is PBST containing 0.05% Tween 20 and 5% skim milk; suggested incubation time is 1 hour at room temperature). Optimal conditions must be determined individually for each application.
Expiration Date	12 months from date of receipt.