

Product Datasheet

BimS/EL/L Antibody [10B12] (orb1272885)

Description

BimS/EL/L Antibody [10B12]

Species/Host

Rat

Reactivity

Canine, Human, Monkey, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, FC, ICC, IHC, IP, WB

Applications
Immunogen

Synthetic peptide corresponding to aa 20-40 of mouse BimL.

Preservatives

Liquid. In PBS containing 0.02% sodium azide and 25% glycerol.

Form/Appearance

Liquid

Concentration

1 mg/ml

Storage

Stable for at least 1 year after receipt when stored at -20°C.

Note

For research use only

Application notes

ELISA: (direct and indirect: 1:2,000-1:5,000). Western Blot: (1:2,000-1:5,000 using ECL; suggested blocking and dilution buffer is PBST containing 0.05% Tween 20 and 5% skim milk; suggested incubation time is 1 hour at room temperature). Optimal conditions must be determined individually for each application.

Isotype

IgG2a, kappa

Clonality

Monoclonal

Dilution Range

ELISA: (direct and indirect: 1:2,000-1:5,000). Western Blot: (1:2,000-1:5,000 using ECL; suggested blocking and dilution buffer is PBST containing 0.05% Tween 20 and 5% skim milk; suggested incubation time is 1 hour at room temperature). Optimal conditions must be determined individually for each application.

Expiration Date

12 months from date of receipt.