

Product Datasheet

TNFSF10 Antibody (orb1272931)

Description	TNFSF10 Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant human TRAIL.
Target	TNFSF10
Preservatives	Liquid. In PBS containing 10% glycerol and 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	1 mg/ml
Storage	Stable for at least 1 year after receipt when stored at -20°C.
Note	For research use only
Application notes	Western Blot: (1:1000).
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P50591
NCBI	4507593
Dilution Range	Western Blot: (1:1000).
Expiration Date	12 months from date of receipt.