



Product Datasheet

TRERF1 antibody (orb127299)

Description	Rabbit Polyclonal to TRERF1.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	C-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q96PN7
Dilution Range	IHC-P: 1:50-1:100, ELISA: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from c...



Immunohistochemical analysis of formalin...