

Product Datasheet

SARS-CoV-2 (COVID-19) Nucleocapsid Antibody [3868] (orb1273019)

Description

SARS-CoV-2 (COVID-19) Nucleocapsid Antibody

Species/Host

Mouse

Reactivity

Virus

Conjugation

Unconjugated

Tested

ELISA, IF

Applications

Immunogen

These newer COVID-19 nucleoprotein Mabs were developed using a recombinant, full-length nucleoprotein and then confirmed that they are reactive with the native virus itself.

Preservatives

Coronavirus Nucleocapsid Antibody is in a phosphate saline buffer (0.01M, pH 7.2) containing 0.1% sodium azide preservative. No stabilizing proteins have been added.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Coronavirus Nucleocapsid Antibody can be stored at -20°C, stable for one year.

Note

For research use only

Application notes

Suggested starting ranges are 1:10-1:50 for IFA and 1:20-1:200 for ELISA. BEST ELISA pairs (capture-conjugate): clones 3865-3868 / 3869-3868

Isotype

Mouse IgG1

Clonality


Monoclonal

Dilution Range


Suggested starting ranges are 1:10-1:50 for IFA and 1:20-1:200 for ELISA. BEST ELISA pairs (capture-conjugate): clones 3865-3868 / 3869-3868

Expiration Date

12 months from date of receipt.



Sandwich
ELISA for
SARS-CoV-
2 (COVID-
19)...



Sandwich
ELISA for
SARS-CoV-
2 (COVID-
19)...