

Product Datasheet

SARS-CoV-2 (COVID-19) Nucleocapsid Antibody [3867] (orb1273020)

Description SARS-CoV-2 (COVID-19) Nucleocapsid Antibody [3867]

Species/Host Mouse

Reactivity Virus

Conjugation Unconjugated

Tested Applications ELISA, IF

Preservatives Coronavirus Nucleocapsid Antibody is in a phosphate saline buffer (0.01M, pH 7.2) containing 0.1% sodium azide preservative. No stabilizing proteins have been added.

Form/Appearance Liquid

Concentration batch dependent

Storage Coronavirus Nucleocapsid Antibody can be stored at -20°C, stable for one year.

Note For research use only


Application notes Suggested starting ranges are 1:10-1:50 for IFA and 1:20-1:200 for ELISA. BEST ELISA pairs (capture-conjugate): clones 3865-3867

Isotype Mouse IgG2b

Clonality Monoclonal

Dilution Range Suggested starting ranges are 1:10-1:50 for IFA and 1:20-1:200 for ELISA. BEST ELISA pairs (capture-conjugate): clones 3865-3867

Expiration Date 12 months from date of receipt.



Sandwich
ELISA for
SARS-CoV-
2 (COVID-
19)...