


Product Datasheet

SARS-CoV-2 (COVID-19) Nucleocapsid Antibody [3866] (orb1273021)


Description

SARS-CoV-2 (COVID-19) Nucleocapsid


COVID-19 (+) COVID-19 (-)



Immunohistochemistry
Validation of SARS-...



Sandwich ELISA for
SARS-CoV-2 (COVID-
19)...



Sandwich ELISA for
SARS-CoV-2 (COVID-
19)...

Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	These newer COVID-19 nucleoprotein Mabs were developed using a recombinant, full-length nucleoprotein and then confirmed that they are reactive with the native virus itself.
Preservatives	Coronavirus Nucleocapsid Antibody is in a phosphate saline buffer (0.01M, pH 7.2) containing 0.1% sodium azide preservative. No stabilizing proteins have been added.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Coronavirus Nucleocapsid Antibody can be stored at -20°C, stable for one year.
Note	For research use only
Application notes	Suggested starting ranges are 1:10-1:50 for IFA and 1:20-1:200 for ELISA. BEST ELISA pairs (capture-conjugate): clones 3865-3866 / 3869-3866
Isotype	Mouse IgG1
Clonality	Monoclonal
Dilution Range	Suggested starting ranges are 1:10-1:50 for IFA and 1:20-1:200 for ELISA. BEST ELISA pairs (capture-conjugate): clones 3865-3866 / 3869-3866
Expiration Date	12 months from date of receipt.