

## Product Datasheet

# SARS-CoV-2 (COVID-19) Nucleocapsid Antibody [3863] (orb1273030)

**Description** SARS-CoV-2 (COVID-19) Nucleocapsid Antibody [3863]

**Species/Host** Mouse

**Reactivity** Virus

**Conjugation** Unconjugated

**Tested Applications** ELISA, IF

**Preservatives** SARS-CoV-2 (COVID-19, 2019-nCoV) Nucleoprotein antibody is in a phosphate saline buffer (0.01M, pH 7.2) containing 0.1% sodium azide preservative. No stabilizing proteins have been added.

**Form/Appearance** Liquid

**Concentration** batch dependent


**Storage** SARS-CoV-2 COVID-19, 2019-nCoV) Nucleoprotein antibody can be stored at -20°C, stable for one year.

**Note** For research use only


**Isotype** Mouse IgG

**Clonality** Monoclonal

**Expiration Date** 12 months from date of receipt.



Western  
Blot  
Validation  
with  
Human  
Recom...



ELISA Test  
Antibodies:  
SARS-CoV-  
2 Nucleo...