

Product Datasheet

SARS-CoV-2 (COVID-19) Nucleocapsid Antibody [3861] (orb1273032)

Description

SARS-CoV-2 (COVID-19) Nucleocapsid Antibody

KDa 1 2

Species/Host

Mouse

Reactivity

Virus

Conjugation

Unconjugated

Tested

Applications

ELISA, IF, WB

Preservatives

SARS-CoV-2 (COVID-19, 2019-nCoV)
Nucleoprotein antibody is in a phosphate saline buffer (0.01M, pH 7.2) containing 0.1% sodium azide preservative. No stabilizing proteins have been added.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

SARS-CoV-2 COVID-19, 2019-nCoV)
Nucleoprotein antibody can be stored at -20°C, stable for one year.

Note

For research use only

Isotype


Mouse IgG1

Clonality


Monoclonal

Expiration Date


12 months from date of receipt.



Western
Blot
Validation
with
Human
Recom...



ELISA Test.
Antibodies:
SARS-CoV-
2 Nucle...



ELISA Test.
Antibodies:
SARS-CoV-
2 Nucle...