

Product Datasheet

SARS-CoV-2 (COVID-19) Nucleocapsid Antibody [3851] (orb1273033)

Description SARS-CoV-2 (COVID-19) Nucleocapsid Antibody [3851]

Species/Host Mouse

Reactivity Virus

Conjugation Unconjugated

Tested Applications ELISA, IF, WB

Preservatives SARS-CoV-2 (COVID-19, 2019-nCoV) Nucleoprotein antibody is in a phosphate saline buffer (0.01M, pH 7.2) containing 0.1% sodium azide preservative. No stabilizing proteins have been added.

Form/Appearance Liquid

Concentration batch dependent

Storage SARS-CoV-2 COVID-19, 2019-nCoV) Nucleoprotein antibody can be stored at -20°C, stable for one year.

Note For research use only


Isotype Mouse IgG2b

Clonality Monoclonal

Expiration Date 12 months from date of receipt.



Western
Blot
Validation
with
Human
Recom...



ELISA Test.
Antibodies:
SARS-CoV-
2 Nucle...