

Product Datasheet

Diphtheria toxin Antibody (FITC) (orb1273422)

Description	Diphtheria toxin Antibody (FITC)
Species/Host	Goat
Reactivity	Bacteria
Conjugation	FITC
Tested Applications	WB
Immunogen	Purified diphtheria toxoid
Preservatives	Diphtheria toxin antibody is in a phosphate saline buffer (0.01M, pH 7.2) containing bovine serum albumin (10 mg/ml) as a stabilizer and 0.1% sodium azide preservative.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Diphtheria toxin antibody can be stored at -20°C, stable for one year.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.