


Product Datasheet

CHAT Antibody (orb1273480)

Description	CHAT Antibody
Species/Host	Goat
Reactivity	Gallus, Goat, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Native choline acetyltransferase purified from human placenta.
Target	CHAT
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Choline Acetyltransferase antibody can be stored at -20°C and is stable at -20°C for at least 1 year.
Note	For research use only
Clonality	Polyclonal
MW	70
Uniprot ID	P28329
NCBI	P28329
Expiration Date	12 months from date of receipt.



Western blot of rat brain lysate showing...