


Product Datasheet

Mitochondria Antibody [113-1] (orb1273646)

Description	Mitochondria Antibody [113-1]
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Human cell homogenate
Form/Appearance	Liquid
Concentration	batch dependent
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	60
Expiration Date	12 months from date of receipt.



Western blot of HeLa lysate showing specific immunolabeling of the ~60k mitochondrial protein.