



Product Datasheet

ZCCHC3 Antibody (orb1273673)

Description	ZCCHC3 Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Ni-NTA purified truncated recombinant human trim5 alpha expressed in E. Coli strain BL21 (DE3).
Target	ZCCHC3
Preservatives	Ascitic fluid containing 0.03% sodium azide.
Concentration	batch dependent
Storage	Trim5α monoclonal antibody can be stored at -20°C, stable for one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	Western Blot:Dilution 1:1,000 - 1:5,000 IHC(P):Dilution 1:100 - 1:400IHC(F):Dilution 1:200 - 1:400Determining optimal working dilutions by titration test.
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q9NUD5
NCBI	Q9NUD5
Dilution Range	Western Blot:Dilution 1:1,000 - 1:5,000 IHC(P):Dilution 1:100 - 1:400IHC(F):Dilution 1:200 - 1:400Determining optimal working dilutions by titration test.
Expiration Date	12 months from date of receipt.



Western blot analysis using Trim5α monoc...



Western blot analysis using Trim5α monoc...