



Product Datasheet

SORL1 Antibody (orb1273685)

Description	SORL1 Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Purified truncated recombinant SORL1-MBP expressed in E. Coli strain BL21 (DE3).
Target	SORL1
Preservatives	Ascitic fluid containing 0.03% sodium azide.
Concentration	batch dependent
Storage	SORL1 monoclonal antibody can be stored at -20°C, stable for one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	Western Blot:1:500 - 1:2,000 IHC(P):1:500 - 1:2,000ELISA:Propose dilution 1:10,000.Determining optimal working dilutions by titration test.
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q92673
NCBI	Q92673
Dilution Range	Western Blot:1:500 - 1:2,000 IHC(P):1:500 - 1:2,000ELISA:Propose dilution 1:10,000.Determining optimal working dilutions by titration test.
Expiration Date	12 months from date of receipt.



Western blot analysis using SORL1 monoclonal antibody.



Immunohistochemical analysis of paraffin-embedded tissue sections.