



Product Datasheet

SNCG Antibody (orb1273686)

Description	SNCG Antibody	
Species/Host	Mouse	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, IHC, WB	
Immunogen	Ni-NTA purified truncated recombinant SNCG expressed in E. Coli strain BL21 (DE3).	
Target	SNCG	Western blot analysis using SNCG monoclo...  Immunohistochemical analysis of paraffin...
Preservatives	Ascitic fluid containing 0.03% sodium azide.	
Concentration	batch dependent	
Storage	SNCG (breast cancer-specific protein 1) monoclonal antibody can be stored at -20°C, stable for one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	Western Blot:1:500 - 1:2,000 IHC(P):1:500 - 1:2,000IHC(F):1:500 - 1:2,000ELISA:Propose dilution 1:10,000.Determining optimal working dilutions by titration test.	
Isotype	IgG1	
Clonality	Monoclonal	
Uniprot ID	O76070	
NCBI	O76070	
Dilution Range	Western Blot:1:500 - 1:2,000 IHC(P):1:500 - 1:2,000IHC(F):1:500 - 1:2,000ELISA:Propose dilution 1:10,000.Determining optimal working dilutions by titration test.	
Expiration Date	12 months from date of receipt.	