

Product Datasheet

RAB25 Antibody (orb1273694)

Description

RAB25 Antibody

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

Ni-NTA purified recombinant Rab25 expressed in E. Coli strain BL21 (DE3).

Target

RAB25

Preservatives

Purified antibody in PBS containing 0.03% sodium azide.

Concentration

batch dependent

Storage

Rab25 monoclonal antibody can be stored at -20°C, stable for one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

Western Blot:1:500 - 1:2,000 IHC(P):1:500 - 1:2,000IHC(F):1:500 - 1:2,000ELISA:Propose dilution 1:10,000.Determining optimal working dilutions by titration test.

Isotype

IgG2a

Clonality

Monoclonal

Uniprot ID

[P57735](#)

NCBI

[P57735](#)


Dilution Range

Western Blot:1:500 - 1:2,000 IHC(P):1:500 - 1:2,000IHC(F):1:500 - 1:2,000ELISA:Propose dilution 1:10,000.Determining optimal working dilutions by titration test.


Expiration Date

12 months from date of receipt.

1 2 3 4



Western blot analysis using Rab25 monoclonal antibody.



Immunohistochemical analysis of paraffin-embedded tissue sections.