



Product Datasheet

MCL1 Antibody (orb1273709)

Description	MCL1 Antibody	
Species/Host	Mouse	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	ELISA, IHC, WB	
Immunogen	Ni-NTA purified recombinant protein MCL-1 expressed in E. Coli strain M15.	
Target	MCL1	Western blot analysis using MCL-1 monocl...  Immunohistochemical analysis of paraffin...
Preservatives	Purified antibody in PBS containing 0.03% sodium azide.	
Concentration	batch dependent	
Storage	MCL-1 monoclonal antibody can be stored at -20°C, stable for one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	Western Blot:Dilution 1:200 - 1:1,000 IHC(P):Dilution 1:200 - 1:1,000 IHC(F):Dilution 1:200 - 1:1,000 ELISA:Propose dilution 1:10,000.Determining optimal working dilutions by titration test.	
Isotype	IgG1	
Clonality	Monoclonal	
Uniprot ID	Q07820	
NCBI	NP_001184249.1	
Dilution Range	Western Blot:Dilution 1:200 - 1:1,000 IHC(P):Dilution 1:200 - 1:1,000 IHC(F):Dilution 1:200 - 1:1,000 ELISA:Propose dilution 1:10,000.Determining optimal working dilutions by titration test.	
Expiration Date	12 months from date of receipt.	