



Product Datasheet

IGF1R Antibody (orb1273723)

Description	IGF1R Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Ni-NTA purified truncated recombinant IGF1R-His expressed in E. Coli strain BL21 (DE3).
Target	IGF1R
Preservatives	Ascitic fluid containing 0.03% sodium azide.
Concentration	batch dependent
Storage	IGF1R monoclonal antibody can be stored at -20°C, stable for one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	Western Blot:1:500 - 1:2,000 IHC(P):1:500 - 1:2,000 ELISA:Propose dilution 1:10,000. Determining optimal working dilutions by titration test.
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	P08069
NCBI	P08069
Dilution Range	Western Blot:1:500 - 1:2,000 IHC(P):1:500 - 1:2,000 ELISA:Propose dilution 1:10,000. Determining optimal working dilutions by titration test.
Expiration Date	12 months from date of receipt.



Western blot analysis using IGF1R monoclonal antibody.



Immunohistochemical analysis of paraffin-embedded tissue sections.