



Product Datasheet

GSK3A Antibody (orb1273733)

Description	GSK3A Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Ni-NTA purified truncated recombinant GSK3 alpha expressed in E. Coli strain BL21 (DE3).
Target	GSK3A
Preservatives	Ascitic fluid containing 0.03% sodium azide.
Concentration	batch dependent
Storage	GSK3 alpha monoclonal antibody can be stored at -20°C, stable for one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	Western Blot:1:500 - 1:2,000 IF(1:200 - 1:1,000 ELISA:Propose dilution 1:10,000.Determining optimal working dilutions by titration test.
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P49840
NCBI	P49840
Dilution Range	Western Blot:1:500 - 1:2,000 IF(1:200 - 1:1,000 ELISA:Propose dilution 1:10,000.Determining optimal working dilutions by titration test.
Expiration Date	12 months from date of receipt.



Western blot analysis using GSK3 monoclo...



Immunofluorescence analysis of HeLa cells...