


Product Datasheet

gfp Antibody (orb1273735)

Description	gfp Antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Ni-NTA purified full-length recombinant GFP-His expressed in E. Coli strain BL21 (DE3).
Target	gfp
Preservatives	Purified antibody in PBS containing 0.03% sodium azide.
Concentration	batch dependent
Storage	GFP monoclonal antibody can be stored at -20°C, stable for one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	Western Blot:1:500 - 1:2,000 ELISA:Propose dilution 1:10,000.Determining optimal working dilutions by titration test.
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	B6UPG7
NCBI	YP_002302326.1
Dilution Range	Western Blot:1:500 - 1:2,000 ELISA:Propose dilution 1:10,000.Determining optimal working dilutions by titration test.
Expiration Date	12 months from date of receipt.



Western blot analysis using GFP monoclon...