

Product Datasheet

CSNK1A1 Antibody (orb1273754)

Description

CSNK1A1 Antibody

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

Ni-NTA purified truncated recombinant CK1 expressed in E. Coli strain BL21 (DE3).

Target

CSNK1A1

Preservatives

Ascitic fluid containing 0.03% sodium azide.

Concentration

batch dependent

Storage

CK1 monoclonal antibody can be stored at -20°C, stable for one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

Western Blot:1:500 - 1:2,000 IHC(P):1:500 - 1:2,000IHC(F):1:500 - 1:2,000ELISA:Propose dilution 1:10,000.Determining optimal working dilutions by titration test.

Isotype

IgG1

Clonality

Monoclonal

Uniprot ID

[P48729](#)

NCBI


[P48729](#)

Dilution Range


Western Blot:1:500 - 1:2,000 IHC(P):1:500 - 1:2,000IHC(F):1:500 - 1:2,000ELISA:Propose dilution 1:10,000.Determining optimal working dilutions by titration test.

Expiration Date

12 months from date of receipt.



Western blot analysis using CK1 monoclon...



Immunohistochemical analysis of paraffin...