



Product Datasheet

PECAM1 Antibody (orb1273760)

Description	PECAM1 Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Ni-NTA purified recombinant CD31 expressed in E. Coli strain BL21 (DE3).
Target	PECAM1
Preservatives	Ascitic fluid containing 0.03% sodium azide.
Concentration	batch dependent
Storage	CD31 monoclonal antibody can be stored at -20°C, stable for one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	Western Blot:1:500 - 1:2,000 IHC(P):1:500 - 1:2,000ELISA:Propose dilution 1:10,000.Determining optimal working dilutions by titration test.
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P16284
NCBI	P16284
Dilution Range	Western Blot:1:500 - 1:2,000 IHC(P):1:500 - 1:2,000ELISA:Propose dilution 1:10,000.Determining optimal working dilutions by titration test.
Expiration Date	12 months from date of receipt.



Western blot analysis using CD31 monoclo...



Immunohistochemical analysis of paraffin...