

Product Datasheet

NPPB Antibody (orb1273772)

Description

NPPB Antibody

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC

Applications

Immunohistochemical
analysis of paraffin...

Immunogen

Synthetic peptide corresponding to aa (Cys-Phe-Gly-Arg-Lys-Met-Asp-Arg-Ile-Ser-Ser-Ser-Ser) of human BNP1

Target

NPPB

Preservatives

Ascitic fluid containing 0.03% sodium azide.

Concentration

batch dependent

Storage

BNP1 monoclonal antibody can be stored at -20°C, stable for one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

IHC(P):1:500 - 1:2,000ELISA:Propose dilution 1:10,000.Determining optimal working dilutions by titration test.

Isotype

IgG1

Clonality

Monoclonal

Uniprot ID

P16860

NCBI

P16860

Dilution Range

IHC(P):1:500 - 1:2,000ELISA:Propose dilution 1:10,000.Determining optimal working dilutions by titration test.

Expiration Date

12 months from date of receipt.