

Product Datasheet

APOM Antibody (orb1273778)

Description

APOM Antibody

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

Ni-NTA purified full-length recombinant ApoM.trx expressed in E. Coli strain BL21 (DE3).

Target

APOM

Preservatives

Purified antibody in PBS containing 0.03% sodium azide.

Concentration

batch dependent

Storage

ApoM monoclonal antibody can be stored at -20°C, stable for one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

Western Blot:1:500 - 1:2,000 IF:1:200 - 1:1,000ELISA:Propose dilution 1:10,000.Determining optimal working dilutions by titration test.

Isotype

IgG1 & IgG2b

Clonality

Monoclonal

Uniprot ID

[O95445](#)

NCBI



[O95445](#)

Dilution Range


Western Blot:1:500 - 1:2,000 IF:1:200 - 1:1,000ELISA:Propose dilution 1:10,000.Determining optimal working dilutions by titration test.

Expiration Date

12 months from date of receipt.



Western blot analysis using ApoM monoclo...



Immunofluorescence analysis of L-02 cell...