

Product Datasheet

DNAJC28 Antibody (orb1273896)

Description	DNAJC28 Antibody
Species/Host	Gallus
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C21orf78.6 protein, 40 kDa
Target	DNAJC28
Form/Appearance	Liquid
Concentration	batch dependent
Storage	2-8°C for one week, -20°C as an undiluted liquid for up to two year.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q9NX36
NCBI	Q9NX36
Expiration Date	12 months from date of receipt.