

Product Datasheet

NT3 Antibody (orb1273910)

Description	NT3 Antibody
Species/Host	Gallus
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	NT3 Mixed Peptides
Target	NT3
Form/Appearance	Liquid
Concentration	batch dependent
Storage	2-8°C for one week, -20°C as an undiluted liquid for up to two year.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.