

Product Datasheet

Phospho Serine Antibody (orb1273921)

Description	Phospho Serine Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IHC, WB
Immunogen	Phosphoserine-KLH and phosvitin mixture
Preservatives	PBS Buffer containing 0.02% sodium azide.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	2-8°C for one week and -20°C as an undiluted liquid for up to two years.
Note	For research use only
Application notes	Phospho-Serine antibody can be used at 1:500 to 1:1,000 for ELISA, 1:200-1:1,000 for Western Blot testing and 1ug/mg sample for IP. Customer using ECL system may need to further dilute the Abs. User needs to try out the titration for optimal concentration.
Clonality	Polyclonal
Dilution Range	Phospho-Serine antibody can be used at 1:500 to 1:1,000 for ELISA, 1:200-1:1,000 for Western Blot testing and 1ug/mg sample for IP. Customer using ECL system may need to further dilute the Abs. User needs to try out the titration for optimal concentration.
Expiration Date	12 months from date of receipt.