

Product Datasheet

SARS-CoV-2 (COVID-19) Spike RBD Monoclonal Antibody [B002] (Biotin) (orb1274007)

Description SARS-CoV-2 (COVID-19) Spike RBD Monoclonal Antibody [B002] (Biotin)

Species/Host Human

Reactivity Virus

Conjugation Biotin

Tested Applications ELISA, NeA, WB

Immunogen S-RBD recombinant protein

Target S

Preservatives PBS with 0.02% sodium azide, 10% glycerol, pH 7.2

Form/Appearance Liquid

Concentration batch dependent

Storage Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.

Note For research use only


Application notes WB: 1:1000~2000

Clonality Monoclonal

NCBI [QHD43416.1](#)

Dilution Range WB: 1:1000~2000

Expiration Date 12 months from date of receipt.



Western Blot Validation with Recombinant...