

Product Datasheet

SARS-CoV-2 (COVID-19) Spike S1 Antibody (IN) (orb1274037)

Description SARS-CoV-2 (COVID-19) Spike S1 Antibody (IN)

Species/Host Rabbit

Reactivity Virus

Conjugation Unconjugated

Tested Applications ELISA, WB

Immunogen Synthetic peptide located within the center of SARS-CoV-2 (COVID-19) Spike S1.

Target S

Preservatives PBS with 0.02% sodium azide, 10% glycerol, pH 7.2

Form/Appearance Liquid

Concentration batch dependent

Storage Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.

Note For research use only

Application notes WB: 1:500-1000


Isotype IgG

Clonality Polyclonal

NCBI [QHD43416.1](#)

Dilution Range WB: 1:500-1000

Expiration Date 12 months from date of receipt.



Western Blot
Validation with
Recombinant...

ELISA experiment

	1	2	
A	1.153	0.135	1:1000
B	0.307	0.129	1:10000
C	0.211	0.118	1:20000
D	0.154	0.1	1:40000
E	0.106	0.143	1:60000
F	0.111	0.141	1:160000
G	0.093	0.106	1:320000
H	0.125	0.101	Blank control

ELISA Test.
Coating
original
concentrati...