

## Product Datasheet

### Orf1ab Antibody (orb1274054)

|                            |                                                                                          |
|----------------------------|------------------------------------------------------------------------------------------|
| <b>Description</b>         | Orf1ab Antibody                                                                          |
| <b>Species/Host</b>        | Rabbit                                                                                   |
| <b>Reactivity</b>          | Virus                                                                                    |
| <b>Conjugation</b>         | Unconjugated                                                                             |
| <b>Tested Applications</b> | ELISA, WB                                                                                |
| <b>Immunogen</b>           | Synthetic peptide located within NSP3-IN1 of SARS-CoV-2 (COVID-19) pp1ab.                |
| <b>Target</b>              | Orf1ab                                                                                   |
| <b>Preservatives</b>       | PBS with 0.02% sodium azide, 10% glycerol, pH 7.2                                        |
| <b>Form/Appearance</b>     | Liquid                                                                                   |
| <b>Concentration</b>       | batch dependent                                                                          |
| <b>Storage</b>             | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles. |
| <b>Note</b>                | For research use only                                                                    |
| <b>Application notes</b>   | WB: 1:500-1000                                                                           |
| <b>Isotype</b>             | IgG                                                                                      |
| <b>Clonality</b>           | Polyclonal                                                                               |
| <b>NCBI</b>                | <a href="#">QHD43415.1</a>                                                               |
| <b>Dilution Range</b>      | WB: 1:500-1000                                                                           |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                          |