



Product Datasheet

CDK1 antibody (orb12741)


Description	Rabbit polyclonal antibody to CDK1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	A synthetic phosphorylated of human CDK1(NP_001777.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	34KDa
Uniprot ID	P06493
Dilution Range	WB: 1:200-1:1000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human lung tissue.



Immunohistochemical staining of human prostate tissue.



Western blot analysis of extracts of various cell lines.