

Product Datasheet

SARS-CoV-2 (COVID-19) ORF9B Antibody (orb1274111)

Description

SARS-CoV-2 (COVID-19) ORF9B Antibody

Species/Host

Sheep

Reactivity

Virus

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

GST-SARS-CoV2 ORF9B

Preservatives

Antibody supplied in PBS.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Aliquot and store at -20°C. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

Detection Antibodies: 1 µg/ml Ab, 5% BSA in TBS-T, 4°C overnight Blocking: 5% milk in TBS-T, 1h room temp Secondary Antibodies: 1:15,000 donkey anti-sheep, 1:15,000 donkey anti-rabbit, in TBS-T, 1h room temp Wash: TBS-T, 4x 10 min

Clonality

Polyclonal

Uniprot ID


P0DTD2

Dilution Range

Detection Antibodies: 1 µg/ml Ab, 5% BSA in TBS-T, 4°C overnight Blocking: 5% milk in TBS-T, 1h room temp Secondary Antibodies: 1:15,000 donkey anti-sheep, 1:15,000 donkey anti-rabbit, in TBS-T, 1h room temp Wash: TBS-T, 4x 10 min

Expiration Date

12 months from date of receipt.



Antibody used at 1 µg/mL in WB.