

Product Datasheet

ORF1ab Antibody (orb1274138)

Description	ORF1ab Antibody
Species/Host	Sheep
Reactivity	Virus
Conjugation	Unconjugated
Immunogen	GST-SARS-CoV2 NSP14
Target	ORF1ab
Preservatives	Antibody supplied in PBS.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Aliquot and store at -20°C. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P0DTC1
Expiration Date	12 months from date of receipt.