

Product Datasheet

SARS-CoV-2 (COVID-19) Nucleocapsid Antibody [6A8] (orb1274167)

Description	SARS-CoV-2 (COVID-19) Nucleocapsid Antibody [6A8]
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant SARS-CoV-2 N protein (His-tag)
Target	N
Preservatives	0.01M PBS(pH7.4) with 0.1% Proclin300
Form/Appearance	Liquid
Concentration	≥1 mg/mL
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	WB(1:500-2000), ELISA(1:5000-1000)
Isotype	IgG2b
Clonality	Monoclonal
Dilution Range	WB(1:500-2000), ELISA(1:5000-1000)
Expiration Date	12 months from date of receipt.